



foodproof® Product List

order no.	real-time PCR Kits	applicable to				
		LightCycler				AriaMx, PikoReal, ABI 7500 and others*
		1.x	2.0	480	96	

foodproof® Foodborne Pathogens

R 302 04	foodproof® <i>Brucella</i> Detection Kit (96 rxn)	-	-	+	+	+
R 302 05	foodproof® <i>Campylobacter</i> Quantification Kit (96 rxn)	-	-	+	+	+
R 310 05	foodproof® <i>Campylobacter</i> Detection Kit (96 rxn)	+	+	+	-	-
R 602 40	foodproof® <i>Clostridium botulinum</i> Detection LyoKit (96 rxn)	-	+	+	+	+
RDK 100 23	foodproof® <i>Clostridium</i> Control Detection Kit (64 rxn)	-	+	+	+	+
R 300 09	foodproof® <i>E. coli</i> and <i>Shigella</i> Detection Kit (96 rxn)	+	+	+	-	-
R 300 10	foodproof® <i>E. coli</i> O157 Detection Kit (96 rxn) (AOAC)	+	+	-	-	-
R 302 10	foodproof® <i>E. coli</i> O157 Detection Kit (96 rxn) (AOAC)	-	-	+	+	+
R 302 15.1	foodproof® <i>Enterobacteriaceae</i> plus <i>Cronobacter</i> Detection Kit (96 rxn) (MicroVal)	-	-	+	+	+
R 310 15.1	foodproof® <i>Enterobacteriaceae</i> plus <i>Cronobacter</i> Detection Kit (96 rxn) (MicroVal)	-	+	+	-	-
R 300 20	foodproof® <i>Listeria</i> Genus Detection Kit (96 rxn)	+	+	-	-	-
R 302 20	foodproof® <i>Listeria</i> Genus Detection Kit (96 rxn)	-	-	+	+	+
R 300 23	foodproof® <i>Listeria monocytogenes</i> Detection Kit (96 rxn) (AOAC, NordVal)	+	+	-	-	-
R 302 23	foodproof® <i>Listeria monocytogenes</i> Detection Kit (96 rxn) (AOAC, NordVal)	-	-	+	+	+
R 602 23	foodproof® <i>Listeria monocytogenes</i> Detection LyoKit (96 rxn)	-	-	+	+	+
R 302 38	foodproof® Norovirus Detection Kit (GI, GII, GIV) (64 rxn) NEW!	-	-	+	+	+
R 310 27	foodproof® <i>Salmonella</i> Detection Kit (96 rxn) (AOAC, MicroVal)	+	+	+	-	-
R 302 27	foodproof® <i>Salmonella</i> Detection Kit (96 rxn) (AOAC, MicroVal, NordVal)	-	-	+	+	+
R 302 27 L	foodproof® <i>Salmonella</i> Detection Kit Large (480 rxn) (AOAC, MicroVal, NordVal)	-	-	+	+	+
R 602 27	foodproof® <i>Salmonella</i> Detection LyoKit (96 rxn) (AOAC)	-	-	+	+	+
R 602 34	foodproof® <i>Salmonella</i> Enteritidis & Typhimurium Detection LyoKit (96 rxn)	-	-	+	+	+
R 602 11	foodproof® STEC Screening LyoKit (96 rxn)	-	-	+	+	+
R 602 12	foodproof® STEC Identification LyoKit (48 rxn) NEW!	-	-	+	+	+
R 602 44	foodproof® <i>Vibrio</i> Detection LyoKit (96 rxn) NEW!	-	-	+	+	+
RDK 100 11	foodproof® <i>Yersinia enterocolitica</i> Detection Kit (64 rxn) NEW!	-	-	+	+	+

* Please request

Note: RDK = Pre-Launch Version

foodproof® is a registered trademark of BIOTECON Diagnostics. The above mentioned real-time PCR instruments are registered trademarks of their respective holders.



foodproof® GMO

RDK 302 39	foodproof® CaMV Detection Kit (96 rxn)	-	-	+	+	+
R 302 29	foodproof® GMO 35S Maize Quantification Kit (2 x 64 rxn)	-	-	+	+	+
R 302 16	foodproof® GMO Bt176 Maize Quantification Kit (2 x 64 rxn)	+	+	+	+	+
R 300 19	foodproof® GMO RR Soya Quantification Kit (2 x 64 rxn)	+	+	+	-	-
R 302 19	foodproof® GMO RR Soya Quantification Kit (2 x 64 rxn)	-	-	+	+	+
R 302 35	foodproof® GMO RR 2 Yield Soya Quantification Kit (2 x 64 rxn)	-	-	+	+	+
R 300 17	foodproof® GMO Screening Kit (35S, NOS) (2 x 64 rxn)	+	+	-	-	-
R 302 17	foodproof® GMO Screening Kit (35S, NOS, bar, FMV) (2 x 64 rxn)	-	-	+	+	+
R 602 17	foodproof® GMO Screening 1 LyoKit (35S, NOS, FMV) (96 rxn) NEW!	-	-	+	+	+
R 602 18	foodproof® GMO Screening 2 LyoKit (5 additional targets) (96 rxn) NEW!	-	-	+	+	+
R 602 21	foodproof® Plant Detection LyoKit (96 rxn) NEW!	-	-	+	+	+

foodproof® Animal Identification

RDK 100 07	foodproof® Donkey Detection Kit (64 rxn)	-	-	+	+	+
R 602 43	foodproof® Porcine Detection LyoKit (96 rxn) NEW!	-	-	+	+	+

foodproof® Spoilage Organisms

R 302 28	foodproof® <i>Alicyclobacillus</i> Detection Kit (96 rxn)	-	-	+	+	+
R 310 02	foodproof® Beer Screening Kit (96 rxn)	+	+	+	-	-
R 302 08	foodproof® <i>Dekkera</i> Quantification Kit (96 rxn)	-	-	+	+	+
R 300 26	foodproof® <i>Saccharomyces cerevisiae</i> var. <i>diastaticus</i> Detection Kit (96 rxn)	+	+	-	-	-
RDK 100 13	foodproof® Yeast and Mold Quantification LyoKit (96 rxn) NEW!	-	-	+	+	+

foodproof® Allergens

R 302 60	foodproof® Celery Detection Kit (64 rxn)	-	-	+	+	+
R 302 61	foodproof® Soya Detection Kit (64 rxn)	-	-	+	+	+
R 302 62	foodproof® Hazelnut Detection Kit (64 rxn)	-	-	+	+	+
R 302 63	foodproof® Peanut Detection Kit (64 rxn)	-	-	+	+	+
R 302 64	foodproof® Gluten Detection Kit (64 rxn) NEW!	-	-	+	+	+

Note: RDK = Pre-Launch Version

foodproof® is a registered trademark of BIOTECON Diagnostics. The above mentioned real-time PCR instruments are registered trademarks of their respective holders.



foodproof® Product List

order no.	DNA Extraction Kits	applicable to
S 400 01	foodproof® ShortPrep I Kit (96 rxn) (AOAC, NordVal)	<i>Salmonella</i>
S 400 02	foodproof® ShortPrep II Kit (96 rxn) (AOAC, NordVal)	<i>Listeria, Campylobacter, E. coli, Shigella</i>
S 400 03	foodproof® ShortPrep III Kit (96 rxn)	Beer spoilage organisms
S 400 04	foodproof® Sample Preparation Kit I (100 rxn)	<i>Enterobacteriaceae, Salmonella, Cronobacter, E. coli, Shigella</i>
S 400 05	foodproof® Sample Preparation Kit II (100 rxn)	<i>Listeria, Campylobacter, Alicyclobacillus</i>
S 400 06.1	foodproof® Sample Preparation Kit III (50 rxn)	GMOs, Allergens and Animal ID
S 400 07	foodproof® StarPrep One Kit (100 rxn) (AOAC, MicroVal)	<i>Enterobacteriaceae, Salmonella, Cronobacter, E. coli, Shigella</i>
S 400 07 L	foodproof® StarPrep One Kit Large (500 rxn) (AOAC, MicroVal)	<i>Enterobacteriaceae, Salmonella, Cronobacter, E. coli, Shigella</i>
S 400 08	foodproof® StarPrep Two Kit (96 rxn)	<i>Listeria, Campylobacter, Alicyclobacillus, Yeast & Mold</i>
S 400 18	foodproof® StarPrep Three Kit (96 rxn)	Beer, <i>Vibrio</i> , STEC
S 400 09	foodproof® StarPrep Four Kit (96 rxn)	Yeast
S 400 11	foodproof® Magnetic Preparation Kit I (96 rxn) (AOAC, MicroVal)	<i>Salmonella, Cronobacter, E. coli, Shigella</i>
S 400 11 L	foodproof® Magnetic Preparation Kit I Large (480 rxn) (AOAC, MicroVal)	<i>Salmonella, Cronobacter, E. coli, Shigella</i>
S 400 12	foodproof® Magnetic Preparation Kit II (96 rxn)	<i>Listeria</i>
S 400 12 L	foodproof® Magnetic Preparation Kit II Large (480 rxn)	<i>Listeria</i>
S 400 13 L	foodproof® Magnetic Preparation Kit III Large (480 rxn)	GMO, Allergens and Animal ID
RDK 400 14	foodproof® Virus Sample Preparation Kit (50 rxn) NEW!	Virus
S 400 15	foodproof® Magnetic Preparation Kit IV (96 rxn) NEW!	<i>Enterobacteriaceae, Salmonella, Cronobacter</i>
S 400 15 L	foodproof® Magnetic Preparation Kit IV (480 rxn) NEW!	<i>Enterobacteriaceae, Salmonella, Cronobacter</i>
S 400 19 L	foodproof® Magnetic Preparation Kit V (480 rxn) NEW!	<i>Salmonella</i> in PPS samples

order no.	Additional Components	applicable to
A 500 02	Reagent D (30 ml/100 rxn) (MicroVal)	foodproof® <i>Enterobacteriaceae</i> plus <i>Cronobacter</i> Detection Kit, foodproof® Yeast and Mold Quantification LyoKit
A 500 02 S	Reagent D (15 ml/100 rxn) (MicroVal)	foodproof® <i>Vibrio</i> Detection LyoKit
A 500 02 L	Reagent D (150 ml/500 rxn) (MicroVal)	foodproof® <i>Enterobacteriaceae</i> plus <i>Cronobacter</i> Detection Kit
A 500 08	Color Compensation Set 1	foodproof® <i>Enterobacteriaceae</i> plus <i>Cronobacter</i> Detection Kit (LightCycler® 2.0, 480)
A 500 10	Color Compensation Set 3	All foodproof® 5' Nuclease kits if used on the LightCycler® 480
A 500 12	Reagent P	foodproof® Magnetic Preparation Kit I
A 500 13	Allergen RM 800	Quantification standard for foodproof® Allergen Detection Kits
A 500 14	Resuspension Reagent	foodproof® StarPrep Two Kit

Note: RDK = Pre-Launch Version

foodproof® is a registered trademark of BIOTECON Diagnostics. The above mentioned real-time PCR instruments are registered trademarks of their respective holders.